

Week: 14


Experiment: 13

“To determine Blood group using the Anti-A, Anti-B, and Anti-D reagents.”

Total Marks (30 Marks):



1) Lab Performance (10 Marks):

2) Observation (5 Marks):



- Anti-A Serum: Dark Blue
- Anti-B Serum: Yellow
- Anti-D (Rh) Serum: More like Translucent solution

3) Picture (2 marks):



4) Result (3 Marks):

ZAIN's Results (Image 1):

- A: Agglutination (clumping) - blue/dark reaction
- D: Agglutination
- B: Agglutination – clumping

Blood Type: AB positive (AB⁺)

EBAD's Results (Image 2):

- A: No clumping- brown/dark reaction
- D: Agglutination - clumping
- B: No clumping - clear/yellow

Blood Type: A positive (A⁺)

5) Labster Simulation (10 Marks):

6) Problems:

Initially we didn't mix it blood anti gens properly, so got wrong results.