

Week: 15

Experiment:14

**“Rapid Detection of Hepatitis-A(HAV) and Dengue Using
Immunochromatographic Test Kits.”**

Total Marks (30 Marks):


1) Lab Performance (10 Marks):

Part A: Detection of Hepatitis-A(HAV)

2) Observations (2 marks):

Lines	Present/Absent
Control line	Present
IgG line	Absent
IgM line	Absent

3) Picture of cassette (1 mark):



4) Results (2 mark):

Negative for Hepatitis (IgG and IgM)

- Control line visible: Test is valid.
- IgM line not visible: No detectable recent infection.
- IgG line not visible: No detectable past infection or long-term antibodies.

Meaning of the result


- You are negative for both IgM and IgG antibodies at this moment.
- This means:
 - No active infection (IgM–)
 - No past resolved infection (IgG–)

Part B: Detection of Dengue

5) Observations (2 marks):

Lines	Present/Absent
Control line	Present
Test line	Absent

6) Picture of cassette (1 mark):



7) Results (2 mark):

Negative for Dengue

- **Control line visible: Test is valid.**
- **NS1 test line not visible: Negative**

So, the test result is considered NEGATIVE.

8) Labster Simulation (10 Marks):

9) Problems:

No such problems faced.